



WITH THE BRAND NEW
DISI evo
DISINTEGRATION TESTER



Tablet & Capsule Testing Equipment To Suit All Your Needs

For Pharmaceutical
Production & Quality Control

www.ischi.com

PRODUCT RANGE



H-Series
Tablet Hardness Tester
H2, H3, H4, H5, H5+

Basic



Hardness



Length/
Diameter



Thickness

Optional



Width



Ext. Scale



Ext.
Thickness
Gauge

[Click here to watch video](#)



P-Series
Semi-Automatic Tablet Testers
P3, P3+, P4, P4+, P5

Basic



Hardness



Thickness



Length/
Diameter

Optional



Width



Weight

[Click here to watch video](#)



Friability + Abrasion Testers
AE-1i with one drum axis

Basic



Friability



Abrasion

Optional

Printer & Balance
for automatic weight
loss calculation

[Click here to watch video](#)



Friability + Abrasion Tester
AE-2i with one drum axis

Basic



Friability



Abrasion

Optional

Printer & Balance
for automatic weight
loss calculation

[Click here to watch video](#)

PRODUCT RANGE



DISI-1M
Manual Disintegration Tester
with 1, 2, 3 or 4 test stations

Basic



Disintegration

[Click here to watch video](#)

NEW



DISI evo

Automatic Disintegration Tester

- Bathless induction heating
- Modular up to 4 stations
- 21 CFR Part II compliant
- Real-time monitoring of the remaining height of the tablets
- Magnetic coupling

Basic



Disintegration

[Click here to watch video](#)



**Automatic Weighing System
for Tablets & Capsules
(available in Lab & IPC versions)**
CIW 7.2/ CIW 7.3

Basic



Weight

Optional



Thickness

[Click here to watch video](#)



CIW6.4
Weight and Length check system
for capsules

Basic



Weight



Length

[Click here to watch video](#)

PRODUCT RANGE



UTS 4.3
Automatic Tablet Testing System
with Touch Software
(available in Lab & IPC Versions)

Basic



Weight



Thickness



Length/
Diameter



Hardness



Length

Optional



Width

[Click here to watch video](#)



UTS-S20
Dust-proof Tablet
Testing System
(OEB level 3)

Basic



Weight



Thickness



Length/
Diameter



Hardness



Width

Optional

[Click here to watch video](#)



UTS X
Automatic test system
with advanced functions
(OEB level 4)

Basic



Weight



Thickness



Length/
Diameter



Hardness

[Click here to watch video](#)



UTS-IP65i
Containment & Wash-in-Place
Tablet Testing System
(OEB level 5)

Basic



Weight



Thickness



Length/
Diameter




Hardness



Protection
IP65

[Click here to watch video](#)

PRODUCT RANGE



Single Punch Tablet Press

This compact and easy to use single punch tablet press machine is able to produce tablet diameters 2–22 mm, upto 50 KN force, producing upto 3,600 tablets per hour, providing unmatched levels of stability, precision and reliability.

[Click here to watch video](#)



Counting + Filling Machine

The professional GMP compliant counting and filling machine for tablets & capsules. Filling and dosing by pure weighing is often too imprecise to achieve consistent results. The benchtop system has a counting range up to 99'999 and a max speed of 57'000/h.

[Click here to watch video](#)

ACCESSORIES



Ph21 Software

For all Testing Devices and Quality Assurance. 21 CFR Part 11 compliance. Enjoy high functionality and full data integrity. Fast and efficient reporting of large amounts of electronic batch data to client servers.

[Click here to go to the website](#)



Dynamic Calibration and Adjustment Device

The dynamic calibration device proves the perfect linearity throughout the whole measurement range and generates a printed calibration report.

The "Dynamic calibration package" comprises:

- Certified external load cell (sensor)
- Software package for the testing device
- Display device

[Click here to go to the website](#)




Mechanical Tablets

The mechanical tablet is a small testing device that allows to conduct quick and easy function control on the hardness measuring station at any time. Each tablet is tested and certified by the manufacturer.

50N, 100N, 150N, 200N

[Click here to go to the website](#)



Turnkey Solution for Tablet Presses

Optimize your tablet production efficiency, by implementing automated tablet testing equipment to reduce costs.

Collect and evaluate much more production data in real time, recognize trends and minimize your risk to manufacture scrap.

[Click here to watch video](#)



Charles Ischi AG
Langfeldstrasse 26, CH 4528 Zuchwil
Switzerland

Tel. +41 (0) 32 621 49 23
E-mail: info@ischi.com
Web: www.ischi.com